



Department of
Primary Industries



Quantitative Analysis of Cooked Rice Grain Texture

Yanco Rice Quality Team


Presented by
Dr Jixun Luo

19th Jul 2017

www.dpi.nsw.gov.au

Introduction

Adapted from DPI website



Rice varieties with different quality and texture are used to make different type of foods

Introduction

Factors affecting cooked rice texture

- Amylose content
- Post-harvest processing
- Cooking methods
 - South and East Asians: rice cooker, particular water ratio
 - Indians: boiling in excess water
 - Americans: use large amount of water then drained


Current methods *(indirectly measures texture)*

- RVA (setback correlated with firmness)
- Gel texture analysis (gel firmness)

Introduction

Texture analysis

- Sensory panel
 - High cost for training and maintaining the panel
 - Not ideal for screening purpose
- Instrumental approach
 - Less cost and less time-consuming
 - Texture analyser
 - Mimic the first bite of a food sample



Aims of study

- Develop a method to directly measure cooked rice texture
- Compare textures of breeding lines and existing varieties
- Determine the contribution of different grain quality factors affecting the texture (*amylose content, gelatinisation temp, RVA parameters, etc.*)

Materials

30 varieties from Leeton farm, C2016

Selected based on

- Different grain dimensions
- Apparent amylose content 7~27 %
- Gelatinisation temp (GT) 65~79 °C
- Optimum cooking time 14~21 min

Method

Cooking method

Excess water method with

- Standard cooking time
- Standard water ratio

Instrumental texture analysis



- Excess water drained after cooking
- Single layer cooked rice grains
- Two cycle compression

Method

Cooking in RVA



- 1g white rice contained in tea ball mesh and rinsed with tap water
- Add water up to 4g in a RVA can
- Seal the can with lid on using the thermo tape
- Load on to RVA machine
- Run the standard cooking profile




- Remove the lid after cooking to release the steam
- Replace with another lid at 1 min to keep it warm for 3 min on the bench

Method

Texture Profile Analysis (TPA) using TXT



- Weigh and prepare 1g cooked rice sample within the mark on a glass plate in 1 min
- Run the standard test profile
 - Two cycle compression, compress to 80% sample height




Technical tips:


Each single grain sits on its side and spread well. Sample height is negatively correlated with firmness for the same sample.

Results - theory


Texture Profile





Results - Firmness




Results - Firmness




Results - Firmness



Results - Stickiness



Results - Stickiness




Amylose
content

7%~9%



17%~21%

25%~27%

Results - Stickiness



Results - Stickiness



Summary

- Used TVT machine to measure texture directly
- Compared results to gel firmness
- Explored texture correlation with OCT, GT, RVA setback parameter
- The outlier varieties in the texture profile of each amylose group were found

Future work

- More quality data for covariant analysis
- Amylose & amylopectin structure analysis by CE/SEC
- Look at the outliers

Acknowledgement

Rice Quality Team:

Dr Rachelle Ward

Ms Margrit Martin

Ms Leanne Johnston

Ms Kylie Elliott

Ms Rachael Wood (PhD candidate)

Dr Laura Pallas

Dr Prakash Oli

Ms Yuki Sims

Dr Mark Talbot

Rice Breeding Team:

Dr Peter Snell

Dr Ben Ovenden

Supervisor:

Dr Sandra McDougall